## LDH MonlabTest®

Piruvato. Cinética UV. DGKC. Líquido.

IVD

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1. Pesce A. Lactate dehydrogenase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1124-117, 438.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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## LDH MonlabTest®

Pyruvate. Kinetic UV. DGKC. Liquid.


### Quantitative determination of lactate dehydrogenase (LDH)

Only for professional in vitro diagnostic use.  
Store at 2-8°C.

#### PRINCIPLE OF THE METHOD

Lactate dehydrogenase (LDH) catalyses the reduction of pyruvate by NADH, according the following reaction:

The rate of decrease in concentration of NADPH, measured photometrically, is proportional to the catalytic concentration of LDH present in the sample<sup>1</sup>.

#### CLINICAL SIGNIFICANCE

Lactate dehydrogenase (LDH) is an enzyme with wide tissue distribution in the body. The higher concentrations of LDH are found in liver, heart, kidney, skeletal muscle and erythrocytes.

Increased levels of the enzyme are found in serum in liver disease, myocardial infarction, renal disease, muscular dystrophy and anemia<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

#### REAGENTS

<b>Reagent 1</b>	Imidazol	65 mmol/L
Buffer	Pyruvate	0,6 mmol/L
<b>Reagent 2</b>	NADH	0,18 mmol/L
Substrate		

#### PRECAUTIONS

R1: H315-Causes skin irritation. H317-May cause an allergic skin reaction. H319-Causes serious eye irritation. Contains: Proclin300.  
Follow the precautionary statements given in MSDS and label of the product.

#### PREPARATION

Working reagent (WR):  
Mix: 4 vol. (R1) Buffer + 1 vol. (R2) Substrate

Stability: 15 days at 2-8°C or 5 days at 15-25°C.

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm <1,00.

#### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C o 37°C (± 0,1°C)
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

#### SAMPLES

Serum<sup>1</sup>. Separated from cells as rapidly as possible. Do not use oxalates as anticoagulants since they inhibit the enzyme.

Do not use haemolysed samples. Stability: 2 days at 2-8°C.

#### PROCEDURE

- Assay conditions:  
Wavelength: .....340 nm  
Cuvette: .....1 cm light path  
Constant temperature: .....25°C /30°C / 37°C
- Adjust the instrument to zero with distilled water or air.
- Pipette into a cuvette:

	25° - 30°C	37°C
WR (mL)	3,0	3,0
Sample (µL)	100	50
- Mix, incubate for 1 minute.
- Read initial absorbance (A) of the sample, start the stopwatch and read absorbance at 1 minute intervals thereafter for 3 minutes.
- Calculate the difference between absorbances and the average absorbance differences per minute (ΔA/min).

#### CALCULATIONS

$$25^{\circ} - 30^{\circ}\text{C} \quad \Delta A/\text{min} \times 4925 = \text{U/L LDH}$$

$$37^{\circ}\text{C} \quad \Delta A/\text{min} \times 9690 = \text{U/L LDH}$$

**Units:** One international unit (IU) is the amount of enzyme that transforms 1 µmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

#### Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,33	1,92
30°C	0,75	1,00	1,43
37°C	0,52	0,70	1,00

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: CONTROL Normal and Pathologic (Ref. MO-165107 and MO-165108).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

#### REFERENCE VALUES<sup>1</sup>

25°C	30°C	37°C
120-240 U/L	160-320 U/L	230-460 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

#### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 3,42 U/L to linearity limit of 1600 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

#### Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	400	785	392	773
SD	3,15	10,97	6,23	9,93
CV (%)	0,79	1,40	1,59	1,28

**Sensitivity:** 1 U/L = 0,00009 ΔA/min.

**Accuracy:** Results obtained using MONLABTEST reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,98382.

Regression equation: y= 0,8988x + 2,583.

The results of the performance characteristics depend on the analyzer used.

#### INTERFERENCES

Haemolysis interferes with the assay.

Some anticoagulants such as oxalates interfere with the reaction<sup>1</sup>.

A list of drugs and other interfering substances with LDH determination has been reported<sup>2,3</sup>.

#### NOTES

**MONLAB has instruction sheets for several automatic analyzers.**

#### BIBLIOGRAPHY

- Pesce A. Lactate dehydrogenase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1124-117, 438.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

#### PACKAGING

Ref.: MO-165092	Ref.: MO-165192
R1: 1 x 60 mL	R1: 1 x 240 mL
R2: 1 x 15 mL	R2: 1 x 60 mL

#### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by