### Calibrador $\beta_2$ -Microglobulina

Monlabtest®

Turbidimetría

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Turbidimetry

#### Quantitative determination of $\beta_2$ -Microglobulin ( $\beta_2$ -m)

**IVD** Only for professional in vitro diagnostic use  
Store at 2-8°C

#### PRODUCT CHARACTERISTICS

The calibrator is a lyophilised serum from human origin used to determinate  $\beta_2$ -m concentration in human samples by turbidimetry or nephelometry methods.

#### REAGENTS

|  |  |
|--|--|
| <b>Calibrator <math>\beta_2</math>-m</b> | Calibrator. Human serum. $\beta_2$ -m concentration is stated on the vial label. |
|--|--|

#### PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

#### CALIBRATION

The concentration value of the Calibrator is calibrated against the 1<sup>st</sup> International  $\beta_2$ -m Standard from WHO.

#### PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Wait for 10 minutes before use.

#### STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 1 month at 2-8°C or 3 months at -20°C.

#### PROCEDURE

To be used in turbidimetric and nephelometric assays.

#### PRESENTATION

|           |          |
|-----------|----------|
| MO-165053 | 1 x 1 mL |
|-----------|----------|