### Calibrador HDL/LDL

Monlabtest®

Suero Humano

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Monlabtest®

Human Serum

#### Quantitative determination of HDL/LDL

**IVD** Only for professional in vitro diagnostic use  
Store at 2-8°C

#### PRODUCT CHARACTERISTICS

HDL/LDL calibrator is a human lyophilised serum. Is intended for determination of HDL or/and LDL concentrations with manual and automated analytical procedures.

#### REAGENTS

Human serum. Biological additives. Bacteriostatics agents. The concentration of the components are lot-specific. See chart below for the values of ingredients.

#### PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

#### CALIBRATION

Information available on request.

#### PREPARATION

- Open one vial carefully, without lost of the lyophilised.
- Reconstitute (→)the lyophilised with 1,0 mL of distilled water.
- Mix thoroughly, avoiding foam forming.
- Bring to room temperature for about 30 min. before use.

Improper handing and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

#### STORAGE AND STABILITY

The Calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial grown.

Store tightly capped when not use.

- After reconstitution is stable for:
  - At 15°C to 25°C 30 hours
  - At 2°C to 8°C 14 days
  - At -25°C to -15°C 3 months

#### PRESENTATION

MO-165111 4 x 1 mL

Component	Method	Value	Units
HDL Cholesterol	Direct-Liquid		mmol/l
LDL Cholesterol	Direct-Liquid		Mmol/l

\*This Control data sheet for lot is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.