**Multicalibrador de Proteínas**  
**Monlabtest®**

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**Multicalibrator Protein**  
**Monlabtest®**

**Quantitative determination of Plasma Proteins**

**IVD** Only for professional in vitro diagnostic use  
Store at 2-8°C

**PRODUCT CHARACTERISTICS**

The Multicalibrator Protein Serum is a liquid serum from human origin containing protein concentrations suitable to calibrate determinations of the constituents listed below by immunoturbidimetry method.

**REAGENTS**

Multicalibrator Protein Serum	Calibrator. Human serum. Sodium azide 0.95 g/L. Concentrations of the constituents are stated on the list below.
-------------------------------	--

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

The concentration values of the Calibrator are traceable to Standard Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM) (all constituents except antithrombin III).

**STORAGE AND STABILITY**

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date. Opened vials, stable for 30 days at 2-8°C.

**PROCEDURE**

To be used in immunoturbidimetric assays.

**PRESENTATION**

MO-165044 1 x 1 mL

PROTEINS	METHOD	CONCENTRATION (mg/dL)
α <sub>1</sub> -antitrypsin	Immunoturbidimetry	365
α <sub>1</sub> -ac. glycoprotein		213
Antithrombin III		78.6
C3		379
C4		71.0
Ceruloplasmin		89.0
Haptoglobin		314
IgA		540
IgG		2780
IgM		238
Prealbumin		76.0
Transferrin		749