Control CK-Nac / CK-MB
MonlabTest®

Suero Control Humano

Suero Control CK/CK-MB

Para uso profesional de diagnóstico *in vitro*.
 Conservar a 2-8°C.

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Human Control Serum

CK/CK-MB Control Serum

Only for professional *in vitro* diagnostic use.
 Store at 2-8°C

PRODUCT CHARACTERISTICS

CK-NAC / CK-MB CONTROL MonlabTest is a human lyophilized control serum.

Is intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Biological additives. Bacteriostatic agents.
 The concentration of the components is lot-specific.
 See chart below for the values of ingredients.

PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

PREPARATION

- Reconstitute (→) with 2.0 mL of distilled water,
 - Mix thoroughly, avoiding foam forming,
 - Bring to room temperature for about 30 minutes before use.
- Improper handling and/or storage can affect results.
 Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:
 - At 15°C to 25°C 24 hours
 - At 2°C to 8°C 3 days
 - At -25°C to -15°C 1 month

Component	Method	T ^a	Value	Range	Units
CK-NAC	DGKC	37 °C			U/L
CK-MB	Immunological Method	37 °C			U/L

PRESENTATION

MO-165110

2 x 2 mL