Multicontrol de Proteínas
Monlabtest®

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Multicontrol Protein
Monlabtest®

Quantitative determination of Plasma Proteins

IVD Only for professional in vitro diagnostic use
Store at 2-8°C

PRODUCT CHARACTERISTICS

The Multicontrol Protein Serum is a liquid serum from human origin containing protein concentrations suitable to evaluate the accuracy and precision of plasma protein determinations by immunoturbidimetry method.

REAGENTS

CONTROL PROTEIN SERUM	Control. Human serum. Sodium azide 0.95 g/L. Protein concentrations are stated on the list below.
------------------------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control are traceable to Standard Reference Material CRM 470/RPHHA (Institute for Reference Materials and Measurements, IRMM) (all constituents except antithrombin III).

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.
Opened vials, stable for 30 days at 2-8°C.

PROCEDURE

To be used in immunoturbidimetric assays.

PRESENTATION

MO-165045 1 x 2 mL

PROTEINS	METHOD	CONCENTRATION	RANGE
α ₁ -antitrypsin	Turbidimetry	65.6	52.5 – 78.7
α ₁ -ac. glycoprotein		50.0	40.0 – 60.0
α ₂ -macroglobulin		102	82 – 122
Antithrombin III		14.1	11.3 – 16.9
C3		78.2	62.6 – 93.8
C4		17.0	13.6 – 20.4
Ceruloplasmin		17.7	14.2 – 21.2
Haptoglobin		83.5	66.8 – 100.2
IgA		155	124 – 186
IgG		655	524 – 786
IgM		79.9	63.9 – 95.9
Prealbumin		16.5	13.2 – 19.8
Transferrin		164	131 – 197