### Control Normal y Patológico

MonlabTest®

IVD

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MonlabTest®

IVD

Serum control multicomponent for clinical chemistry assays

Only for professional *in vitro* diagnostic use.  
Store at 2-8°C.

#### PRODUCT CHARACTERISTICS

Normal or Pathologic Control MonlabTest is a human lyophilised serum. With most constituent concentrations and activities in the normal or pathologic range.

It's intended for control of accuracy for use with manual and automated analytical procedures.

#### REAGENTS

Human serum. Biological additives. Bacteriostatic agents. The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

#### PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

#### PREPARATION

- Reconstitute (→) with 5.0 mL of distilled water,
  - Mix thoroughly, avoiding foam forming,
  - Bring to room temperature for about 30 minutes before use.
- Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

#### STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:
  - At 15°C to 25°C ..... 12 hours
  - At 2°C to 8°C ..... 5 days
  - At -25°C to -15°C ..... 1 month
- Bilirubin (stored protected from light):
  - At 15°C to 25°C ..... 8 - 4 hours (total- direct)
  - At 2°C to 8°C ..... 24 - 8 hours (total- direct)
  - At -25°C to -15°C ..... 2 weeks
- Acid phosphatase:
  - At 15°C to 25°C ..... 4 hours
  - At 2°C to 8°C ..... 1 day
  - At -25°C to -15°C ..... 2 weeks

#### PRESENTATION

Normal Control: MO-165107	4 x 5 mL
Pathological Control: MO-165108	4 x 5 mL