### Control Lipídico Monlabtest®

Suero Control Humano

#### Control Lipídico

**IVD** Para uso profesional de diagnóstico in vitro  
Conservar a 2-8°C

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Human Control Serum

#### Lipidic Control

**IVD** Only for professional in vitro diagnostic use  
Store at 2-8°C

#### PRODUCT CHARACTERISTICS

Lipidic Control Monlabtest is a human lyophilised control serum. Is intended for control of accuracy for use with manual and automated analytical procedures.

#### REAGENTS

Human serum. Biological additives. Bacteriostatics agents. The concentration of the components are lot-specific. See chart below for the values of ingredients.

#### PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

#### CALIBRATION

Information available on request.

#### PREPARATION

- Reconstitute (→) with 1.0 mL of distilled water.
- Mix thoroughly, avoiding foam forming.
- Bring to room temperature for about 30 min. before use.

Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

#### STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:
  - At 15°C to 25°C 30 hours
  - At 2°C to 8°C 14 days
  - At -25°C to -15°C 3 months

#### PRESENTATION

MO-165112 1 x 1 mL

Component	Method	Value	Range	Units
<b>Cholesterol</b>	CHOD-PAP			mmol/l
<b>HDL cholesterol</b>	Direct- Liquid			mmol/l
<b>LDL cholesterol</b>	Direct- Liquid			mmol/l
<b>Phosfolipids</b>	CHOLIN-PAP			mmol/l
<b>Triglycerides</b>	GPO-PAP			mmol/l

\*This Control data sheet for lot is applicable to sublotes. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.