## RPR MonlabTest®

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1. George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
3. Sandra Larsen et al. A manual of Test for Syphilis American Public Health.
4. Joseph Earle Moore et al. Gastrointestinal Hemorrhage 1952; 150(5): 467-473.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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Ref. MO-165020	500 test	Ref. MO-165021	150 test
2 x 5 mL RPR-carbón		3 mL RPR-carbón	
1 mL Control +		1 mL Control +	
1 mL Control -		1 mL Control -	
63 x 8 portas desechables		21 x 8 portas desechables	
Vial y aguja dispensadora		Vial y aguja dispensadora	

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## RPR MonlabTest®

Qualitative determination of plasma reagins

Only for professional *in vitro* diagnostic use. Store at 2-8°C.

### PRINCIPLE OF THE METHOD

The RPR MonlabTest is a non-treponemal slide agglutination test for the qualitative and semi-quantitative detection of plasma reagins in human serum. Carbon particles coated with a lipid complex are agglutinated when mixed with samples containing reagins of patient affected by syphilis.

### CLINICAL SIGNIFICANCE

Reagins are a group of antibodies against some components of the damage tissues from patients infected by *Treponema pallidum*, the agent which causes the syphilis. This microorganism produces some damage to the liver and heart, releasing some tissue fragments. Immunological patient system react producing reagins, antibodies against these fragments. The assay is useful to follow the antibiotic therapy answer.

### REAGENTS

<b>RPR-carbon</b>	Carbon particles coated with a lipid complex, cardiolipin, lecithin and cholesterol in phosphate buffer 20 mmol/L. Preservative. pH 7.0.
<b>Control +</b> Red cap	Artificial serum with reagin titer $\geq 1/4$ .
<b>Control -</b> Blue cap	Animal serum. Preservative

### PRECAUTIONS

Control + / - : H317-May cause an allergic skin reaction. Contain 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product.

### CALIBRATION

The sensitivity is calibrated against the International Reference WHO (1<sup>st</sup> International Standard for Syphilitic Human Serum, ref. 05/132).

### PREPARATION

**RPR-carbon:** Homogenize the reagent before use. Place the needle to the plastic dispenser vial, open the RPR-carbon vial and aspirate the required amount of reagent. Once the test is finished, return the reagent to the original vial and rinse the needle and dispenser with distilled water.

### STORAGE AND STABILITY

All the kit components will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Mix reagents gently before use.  
**Do not freeze:** frozen reagents could change the functionality of the test.  
**Reagents deterioration:** Presence of particles and turbidity.

### ADDITIONAL EQUIPMENT

- Mechanical rotator with adjustable speed at 80-100 r.p.m.
- Vortex mixer.
- Pipettes 50  $\mu$ L.

### SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

### PROCEDURE

#### Qualitative method

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50  $\mu$ L of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
3. Homogenize the reagent before using. Invert the dispenser vial and press lightly to remove air bubbles.
4. Place the dispenser vial together with the needle in a vertical position and perpendicular to the slide and add one drop (20  $\mu$ L) of reagent together with each of the samples and controls.
5. Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
6. Place the slide on a mechanical rotator at 80-100 r.p.m. for 8 min<sup>(Note 1)</sup>. False positive results could appear if the test is read later than 8 minutes.

#### Semi-quantitative method

1. Make serial two fold dilutions of the sample in 9 g/L saline solution.
2. Proceed for each dilution as in the qualitative method.

### READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide test from the rotator. Rotate the slide twice by hand before reading.

#### Interpretation

Agglutination	Reading	Report
Medium or large clumps	R	Reactive
Small clumps	W	Weakly reactive
No clumping or very slight "roughness"	N	Non Reactive

The titer, in the semi-quantitative method, is defined as the highest dilution showing a positive result.

### QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of procedure, as well as a comparative pattern for a better result interpretation. All result different from the negative control result, will be considered as a positive.

### PERFORMANCE CHARACTERISTICS

1. **Analytical sensitivity:** Accurate titer determination of the Reference Material, under the described assay conditions (see calibration).
2. **Prozone effect:** No prozone effect was detected up to titers  $\geq 1/128$ .
3. **Diagnostic sensitivity:** 100%.
4. **Diagnostic specificity:** 100 %.

### INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (10 g/L) and lipids (10 g/L), do not interfere. Rheumatoid factors (300 IU/mL), interfere. Other substances may interfere<sup>5</sup>.

### NOTES

1. During the 8 minutes of reaction time do not expose the slide to a source of heat or intense light in order to reduce evaporation. Such evaporation could cause a false agglutination and therefore false positive results.

### LIMITATIONS OF THE PROCEDURE

- RPR MonlabTest is non-specific for syphilis. All Reactive samples should be retested with treponemic methods such as TPHA and FTA-Abs to confirm the results.
- A Non Reactive result by itself does not exclude a diagnosis of syphilis. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.

### BIBLIOGRAPHY

1. George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
3. Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
4. Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952; 150(5): 467-473.
5. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCPress, 1995.

### PACKAGING

Ref.: MO-165020 RPR 500 tests	Ref.: MO-165021 RPR 150 tests
2 x 5 mL RPR-carbon	3 mL RPR-carbon
1 mL Control +	1 mL Control +
1 mL Control -	1 mL Control -
63 x 8 disposables slides	21 x 8 disposables slides
Dispensing vial and needle	Dispensing vial and needle

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For in vitro diagnostic use only
	Contains sufficient for <n> tests		Consult instructions for use
	Catalogue Code		Keep dry
	Lot Number		Temperature limitation
	Don't re-use		Use by