SARS-CoV-2 MonlabTest®

MO-804032 20 TESTS

One Step to detect SARS-CoV-2 Antigen

IVD

A rapid one step test for the qualitative detection of nucleoprotein antigen of SARS-CoV-2 from human nasopharyngeal swab or human nasal samples specimens.

For professional *in vitro* diagnostic use only.

INTENDED USE

The SARS-CoV-2 MonlabTest® is a rapid chromatographic immunoassay for the qualitative detection of Coronavirus (SARS-CoV-2) antigens in human nasopharyngeal swab or human nasal samples specimens to aid in the diagnosis of Coronavirus (SARS-CoV-2) infection.

SYNTHESIS

Recently, SARS-CoV-2 (Severe Acute Respiratory Syndrome Corona Virus 2) emerged as a new specie able to infect humans and currently causing the COVID-19 (Coronavirus disease-19) pandemic, with an overwhelming hospitalization of infected patients.

Clinically, patients with SAR-CoV-2 infection tend to suffer symptoms such as olfactory and gustatory dysfunction, fever, dry cough, anosmia, fatigue, dyspnoea, headache, diarrhoea and sore throat, followed by vascular and systemic complications, such as leukocyte infiltration of the lungs. COVID-19 commonly results in pneumonia, which can evolve to into acute respiratory distress syndrome, leading to respiratory or multiorgan failure.

Coronaviruses (CoVs) belongs to the order of Nidovirales, identified by its envelope characteristics and positive-sense RNA as genetic material. CoVs consist of four structural proteins: spike (S), membrane (M), envelope (E) and nucleocapsid (N) proteins.

The M and E proteins are essential for virus assembly, while the S protein, on the surface of the viral particles, is crucial for affinity and attachment to host cells. The S protein, responsible for viral entry, determines the host tropism and virus transmission. The N protein, the main structural protein of SARS-CoV-2 is responsible for the transcription and replication of viral RNA, the packaging of the encapsulated genome into virions and interactions with the cell cycle of host cells. In addition, the N protein, which has a substantial immunogenic ability, is abundantly expressed during viral infection.

It is believed that the high transmissibility of COVID-19 is related to its high viral loads in the upper respiratory tract and the fact that many individuals remain asymptomatic, shedding and transmitting the virus. The COVID-19 virus transmission occurs via droplets of saliva or discharge from the nose of infected individuals, therefore following proper hygiene practices when coughing and sneezing are keys to deplete transmission.

PRINCIPLE

The SARS-CoV-2 MonlabTest® is a qualitative lateral flow immunoassay for the detection of Coronavirus (SARS-CoV-2) antigen in human nasopharyngeal samples or human nasal samples. The membrane is pre-coated with monoclonal antibodies against SARS-CoV-2 antigens on the test line region. During testing, the sample reacts with the particle coated with anti-SARS-CoV-2 antibodies which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. In the case of a positive result the specific antibodies present on the membrane will react with the mixture conjugate and generate a coloured line. A green coloured line always appears in the control line and serves as verification that sufficient volume was added that proper flow was obtained and as an internal control for the reagents.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pouch until use.
- Do not use the test if pouch is damaged.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- Clean up spills thoroughly using an appropriate disinfectant.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. Please, take note that the COVID-19 infection is very dangerous, extremely the precautions.
- The test should be discarded in a proper biohazard container after testing.
- The test must be carried out within 2 hours of opening the sealed bag.
- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal or nasal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- The presence of yellow lines in the results window (control and test line zone) that are visible before using the test are completely normal. That does not mean failure on test functionality.
- Positive results should be processed following local laws and regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

MATERIALS PROVIDED

- 20 Tests
- Instructions for use
- 20 Sample diluent vials (sample and positive control diluent)
- 20 Sterile Swabs
- 1 SARS-CoV-2 control +
- 1 Laboratory rack

MATERIALS REQUIRED BUT NO PROVIDED

- Disposable gloves (other personal protective equipment that will be considered necessary)
- Timer
- Vortex (optional)

SPECIMEN COLLECTION AND PREPARATION

Nasopharyngeal swab method:

- Bend shaft to follow curve of nasopharynx.
 - Insert swab through nostril to posterior nasopharynx.
 - Rotate swab a few times to obtain infected cells.
 - Considerer repeating procedure using other nostril only if the protocol of the professional taking the sample required do it.
- Evaluate the specimen immediately (testing sensitivity decrease over time).

Follow immediately with the point test procedure.

Nasopharyngeal samples previously extracted by transport media:

Even while using direct nasopharyngeal swab is the preferred protocol, the device can be also used with samples from nasopharyngeal swabs previously diluted in transport media such as: VTm, UTM or Saline Buffer.

- Use the minimum volume of transport media in order to avoid loss in sensitivity. Example: 1.0 mL. Using higher volumes will affect sensitivity of the system.
- Dilute the extracted sample 1:1 in Reagent (sample diluent) provided. (0.5mL/vial).

Follow immediately with the point test procedure.

Nasal swab method:

- Insert swab through nostril to nasopharynx (approx. 2cm).
- Rotate swab a few times to obtain infected cells.
- Remove the swab from the nostril carefully and repeat the procedure with the other nostril.

Evaluate the specimen immediately (testing sensitivity decrease over time).

Follow immediately with the point test procedure.

Samples should be processed as soon as possible after collection. Follow proper infection control practices.

If test samples after collection is not possible, samples could be stored in refrigerator, 2°-8°C (36°-46.4°F) for a maximum of 8 hours prior to testing. If the samples (only for nasopharyngeal samples) are preserved in validated transport media (VTM, UTM o Saline buffer) could be preserved on it until 6 hours at room temperature or in the refrigerator (2-8°C).

PROCEDURES

Allow the tests, samples, buffer and control to reach room temperature (15-30°C/59-86°F) prior to testing. Do not open pouches until ready to perform the assay.

To process the collected nasopharyngeal swabs/nasal samples (see illustration 1):


Put the tubes in the laboratory rack (identify properly samples from the different patients). Use a separate sample diluent vial for each sample (swab). Open the blue cap and immediately put into the nasopharyngeal or nasal swab, mix for 15 seconds and extract as much liquid possible from the swab. Discard the swab. Close the sample diluent vial (blue cap), be sure that both caps (blue and translucent caps) are properly closed. Samples can be shake using a vortex (optional).

Remove the SARS-CoV-2 MonlabTest® from its sealed pouch and use it as soon as possible. Use a separate device for each sample. Open the translucent cap and dispense 3 drops into the specimen well (S). Start the timer. Read the result at **10 minutes** after dispensing the sample.

The Diluent has SARS-CoV-2 inactivation capability: 99.68% after 1 minute and 99.98% after 10 minutes.


Illustration 1 Nasopharyngeal swab/Nasal swab

Open the blue cap (sample diluent vial)
Put the swab into the sample diluent vial,
Extraction of the sample 15 seconds, then discard the swab and close the blue cap.
Open the translucent cap and add **3 drops** of the mixture into the card (specimen well (S)).

INTERPRETATION OF RESULTS

Illustration 2

POSITIVE: Two lines appear across the central window, a **red** test line marked with the letter T and a **green** control line marked with the letter C.

NEGATIVE: Only one **green** line appears across the control line region marked with the letter C (control line).

INVALID: Total absence of the **green** control coloured line regardless the appearance or not of the **red** test line. Note: Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor. See illustration 2.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured test line in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test:

- A green line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

External Quality Control


Each kit contains a positive control material. Use the control swab to check that the extraction reagents and the test are working properly. Also use the controls to test that you are able to correctly perform the test procedure.

Quality Control Procedure:

SARS-CoV-2 Positive control: Remove the SARS-CoV-2 positive control from its sealed pouch. Open the blue cap and immediately put the SARS-CoV-2 positive control swab into the vial, mix 15 seconds and extract as much liquid possible from the swab. Discard the swab. Close the sample diluent vial (blue cap), be sure that both caps (blue and translucent caps) are properly closed. Remove the test from its sealed pouch, open the translucent cap and dispense 3 drops of the positive control liquid into the specimen well (S).
Result: SARS-CoV-2 positive (see interpretation of results).

Illustration 3. Positive Control

Open the blue cap (sample diluent vial)
Put the **SARS-CoV-2 control +** into the sample diluent vial,
Extraction of the control+ 15 seconds, then discard the swab and close the blue cap.
Open the translucent cap and add **3 drops** of the mixture into the card (specimen well (S)).

LIMITATIONS

1. SARS-CoV-2 MonlabTest® will only indicate the presence of Coronavirus (SARS-CoV-2) in the specimen (qualitative detection) and should be used for the detection of Coronavirus (SARS-CoV-2) antigens in nasopharyngeal specimens and/or nasal specimens (from swab). Neither the quantitative value nor the rate of increase in Coronavirus (SARS-CoV-2) antigens concentration can be determined by this test.
2. Positive results do not rule out co-infections with other pathogens.
3. VTM, UTM and Saline Buffer are the transport media validated for be used with this test (SARS-CoV-2 MonlabTest®). Always, following proportion 1:1 (transport media and sample diluent provided with the test). When using transport media the sensitivity of the test can be reduce due to excessive dilution of sample. This is not recommended for nasal samples. The preference is to use the sample immediately after taking it.
4. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of Coronavirus (SARS-CoV-2) infection.
5. This test provides a presumptive diagnosis of Coronavirus (SARS-CoV-2) respiratory infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

EXPECTED VALUES

On March 11, 2020 the World Health Organization (WHO) declared COVID-19 a pandemic.

COVID-19 has led to an emerging, rapidly evolving situation: reaching more than 200 countries, areas, or territories, to date, and causing thousands of deaths, with an overall mortality around 3%.

Emerging data from patients afflicted by COVID-19 indicate that genetic immunologic, metabolic and environmental factors are involved in the pathogenesis of COVID-19. It was found that certain infected people were asymptomatic. These individuals are an important reservoir for spread of the infection.

PERFORMANCE CHARACTERISTICS

Detection limit

The detection limit value is: 1.25×10^2 PFU/mL of nCoV-2019 D614G(S).

Sensitivity and Specificity

It was performed a multi-center evaluation with 556 nasopharyngeal samples from people suspected of infection by SARS-CoV-2.

The evaluation was performed using SARS-CoV-2 MonlabTest® versus qPCR technique.

IC test: SARS-CoV-2 MonlabTest® (Nasopharyngeal sample)	qPCR technique (nasopharyngeal)			Total
	+	-		
+	93	1		94
-	7	455		462
Total	100	456		556

SARS-CoV-2 MonlabTest® vs qPCR technique		95% CI (Confidence interval)	
Sensitivity	93.0%	86.1-97.1%	
Specificity	99.8%	98.8-100.0%	
PPV	98.9%	94.2-100.0%	
NPV	98.5%	96.9-99.4%	

This multi-center evaluation (nasopharyngeal samples) with positive samples Ct<28 (*) showed following results:

Sensitivity: 95.1 % (95% confidence interval: 88.0-98.7%)

Specificity: 99.8 % (95% confidence interval: 98.8-100.0%)

(*) Recommendations of the WHO: "Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)" were considered. The sensitivity of the test was calculated with nasopharyngeal samples with high viral load in the range of Ag-RDT detection. High viral loads are expected in early symptomatic phases of the illness (first 5-7 days of illness).

It was performed an evaluation with 990 nasal samples from people suspected of infection by SARS-CoV-2. The evaluation was performed using SARS-CoV-2 MonlabTest® versus qPCR Technique.

IC test: SARS-CoV-2 MonlabTest® (Nasal sample)	qPCR technique (nasopharyngeal)			Total
	+	-		
+	129	7		136
-	27	827		854
Total	156	834		990

SARS-CoV-2 MonlabTest® vs qPCR technique		95% CI (Confidence interval)	
Sensitivity	82.7%	75.8-88.3%	
Specificity	99.2%	98.3-99.7%	
PPV	94.9%	89.7-97.9%	
NPV	96.8%	95.4-97.9%	

This evaluation (nasal samples) with positive samples Ct<28 (*) showed following results:

Sensitivity: 96.9 % (95% confidence interval: 91.2-99.4%)

Specificity: 98.5 % (95% confidence interval: 97.7-99.1%)

Hook Effect

SARS-CoV-2 MonlabTest® does not show hook effect at the concentration of SARS-CoV-2 protein tested. (202500.0 ng/mL)

Cross-Reactivity

It was performed an evaluation to determine the cross reactivity of SARS-CoV-2 MonlabTest®. There is not cross reactivity with

common respiratory pathogens, other organisms and other substances that could cause infections:

Adenovirus	Coronavirus strain NL63	Legionella	Respiratory syncytial virus
Astrovirus	Coronavirus strain OC43	Listeria monocytogenes	Salmonella enteritidis/paratyphi A/typhi/typhimurium
Bocavirus	Entamoeba histolytica	Mycobacterium tuberculosis	Shigella boydii/dysenteriae/flexneri/sonnei
Bordetella pertussis	Enterovirus	Mycoplasma pneumoniae	Staphylococcus aureus
Calprotectin (human)	Escherichia coli O157:H7	Mycoplasma tuberculosis	Staphylococcus epidermidis
Campylobacter jejuni	Haemophilus influenzae	Pneumocystis jirovecii	Streptococcus pneumoniae
Chlamydia pneumoniae	Giardia	Metapneumovirus human (hMPV)	Streptococcus pneumoniae
Clostridium difficile GDH/Toxin A /Toxin B	Helicobacter pylori	Parainfluenza virus	Streptococcus pyogenes
Cryptosporidium	Hemoglobin (human/bovine/pig)	Norovirus GI/Norovirus GII	Transferrin (human)
Coronavirus strain 229E	Lactoferrin (human)	Rhinovirus	Yersinia enterocolitica O:3/O:9
Coronavirus strain HKU1	Influenza A/Influenza B	Rotavirus	

SARS-CoV-2 MonlabTest® has some cross reaction with SARS and very low cross reactions with MERS.

Interferences

It was performed an evaluation to determine the possible interferences of SARS-CoV-2 MonlabTest®. There are not interferences with the substances tested:

Metronidazole	Loratadine	Loperamide hydrochloride (Fortasec)	Phenoxymethylpenicillin potassium
Ampicillin	Dexchlorpheniramine (Polaramine)	Heparin (Hitor)	Ambroxol hydrochloride (Mucosan)
Oseltamivir	Ebastine (Ebastel)	Almagato (Almax)	Macrogol 3350 (Movicol)
Amantadine	Acetyl Salicylic (Adiro)	Fosfomicin (Monurol)	Lysine Carbocysteinate (Pectox)
Ribavirin	Ibuprofen (Espidifen)	Acetylcysteine (Flumucil)	Hydroxyzine dihydrochloride
Codeine (Toseina)	Paracetamol (Dolocetil)	Dexketoprofen trometamol (Enantyum)	Lorazepam
Benzocaine (Angleptol)	Metamizole (Nolotil)	Levofloxacin	Amoxicillin
Cloperastine (Flutox)	Prednisone	Ciprofloxacin	Mercaptopurine
Carbocisteine (Iniston mucolítica)	Omeprazole	Sore Throat Phenol spray	Mupirocin
Naso GEL	ZICAM	Homeopathic	Tobramycin
CVS Nasal Spray (Cromolyn)	Afrin (Oxymetazoline)	Biotine	Chloraseptic (Menthyl/Benzocaine)
CVS Nasal Drops (Phenylephrine)	Rifampicin (Rifaldin)	Fluticasone Propionate	Adenovirus
STREP A	RSV	Influenza A	Influenza B

REFERENCES

- Palestino, G.; García-Silva, I.; González-Ortega, O. and Rosales-Mendoza, S. **Can nanotechnology help in the fight against COVID-19? Expert Review of Anti-infective Therapy.** <https://doi.org/10.1080/14787210.2020.1776115>, Jun. 2020.
- Abiodun Iwalokun, B.; Olalekan, A.; Adenipekun, E.; Ojo, O.; Senapon Iwalokun, O.; Mutiu, B.; Orija, O.; Adegbola, R.; Salako, B. and Akinloye, O. **Improving the Understanding of Immunopathogenesis of Lymphopenia as a Correlate of SARS-COV-2 Infection Risk and Disease Progression in African Patients: UGLY SARS-COV-2 Study Protocol.** Journal of Medical Internet Research. <https://preprints.jmir.org/preprint/21242>, Jun. 2020.
- Giménez, E.; Albert, E.; Torres, I.; Remigia, M. J.; Alcaraz, M. J.; Galindo, M. J.; Blasco, M. L.; Solano, C.; Fomer, M. J.; Redón, J.; Signes-Costa, J. and Navarro, D. **SARS-CoV-2-reactive interferon-γ-producing CD8+T cells in patients hospitalized with Coronavirus Disease 2019.** Journal of Medical Virology. <https://doi.org/10.1002/jmv.26213>, Jun. 2020.
- World Health Organization, Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays. Interim guidance. 11 September.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For in vitro diagnostic use only
	Authorized representative		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by
	Diluent (sample diluent)		Do not re-use

RECOMENDATIONS

Recommendations **World Health Organization**: Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays- 11 September 2020.

1. To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within **the first 5-7 days following the onset of symptoms.**

2. Where possible, all positive samples giving positive results should be transported to laboratories with NAAT capability for confirmatory testing.
3. This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT-confirmed COVID-19 outbreaks.
4. To monitor trends in disease incidence in communities.
5. For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
6. A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.
7. Even SARS-CoV-2 MonlabTest® test was not validated using samples from asymptomatic contacts of cases, asymptomatic cases have been demonstrated to have viral loads similar to symptomatic cases, so SARS-CoV-2 MonlabTest® could detect as positive.

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